

Table S1 Equipment details and scan parameters

Center	Equipment	Tube voltage (KV)	Tube current (mA)	Matrix size	Slice thickness(mm)	Slice Spacing(mm)
1	LightSpeed VCT (GE)	120	220	512×512	5	5
2	Iqon Spectral CT (Philips)	120	220	512×512	5	5
2	LightSpeed VCT (GE)	120	220	512×512	5	5
3	Revolution Maxima (GE)	80	320	512×512	5	5
4	SOMATOM definition Flash (SIEMENS)	120	390	512×512	5	5

Center 1, Affiliated Hospital of North Sichuan Medical College, Nanchong, China. Center 2, Chongqing General Hospital, Chongqing, China. Center 3, Chongqing Wulong People's Hospital, Chongqing, China. Center 4, Chongqing Jiangjin Second People's Hospital, Chongqing, China.




Figure S1 Importance of 5 radiomics features of the IL (A), AH (B), and WB (C). IL, infarct lesion; AH, affected hemisphere; WB, whole brain.