

Appendix 1

CT angiography parameters

Manufacturer	Matrix Size	FOV (mm)	Calculated In-Plane Resolution (mm)
Siemens	512×512	332	0.65×0.65 (332/512)
TOSHIBA 320	512×512	227.5	0.44×0.44 (227.5/512)
GE	512×512	365	0.71×0.71 (365/512)

MRI parameters

Sequence	Matrix Size	FOV (mm)	Calculated In-Plane Resolution (mm)
T2-fat suppression	480×480	350	0.73×0.73 (350/480)
Late enhancement	1024×1024	400	0.39×0.39 (400/1024)