

## Supplementary

**Table S1** Magnetic resonance imaging protocol

Sequence	TR (ms)	TE (ms)	FOV (mm <sup>2</sup> )	ST (mm)	Slice gap (mm)
T1WI	678	15	220×209	5	1
T2WI	5,728	110	220×187	5	1
T2 FLAIR	9,480	120	220×188	5	1

T1WI, T1-weighted image; T2WI, T2-weighted image; T2 FLAIR, T2-weighted fluid attenuated inversion recovery; TR, repetition time; TE, echo time; FOV, field of view; ST, slice thickness.